

# HEINE® MINI-C CLIP LAMP

## Compact, long-life pocket diagnostic light

Sturdy, long-life pocket diagnostic light in a compact, modern design with chrome head. Bright, concentrated light thanks to unique XHL Xenon Halogen lens-bulb. mini-c Clip Lamp only fits mini-c handles. The mini-c handle is not compatible with mini 3000 instruments. Handle available in black or blue. Clip Lamp with Ear Light for otoplasty. The Ear Light is used to ensure that hearing aids fit correctly. Compact, oval, attractive light. Lies ergonomically in the hand. High-quality handle: Chrome-finish upper section / refined plastic. Shockproof, sturdy, non-slip. Attachment clip with integrated on / off switch. Secure. Switches off automatically when replaced in the pocket. 20.000 switch cycles guaranteed. Replaceable batteries. Size AAA.



Established 1922

### Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995

