

HEINE® EN 200 BP Automatic Digital Sphygmomanometer

Reliable measurement values due to outstanding measurement accuracy and intuitive operation.

NEW!



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Graphic display with 16 mio. colours and touch screen.

Digital, fully automatic, with intuitive menu navigation. Easy to use, optimises the workflow and saves time.

Capable of memorising recent blood pressure and pulse readings. 500 memory spaces per measurement mode.

Reliably detects **cardiac dysrhythmia and motion artifacts.**



Outstanding measurement accuracy with precise results (+/-3 mmHg, pulse 30–220/min). Meets the high scientific requirements of the European Society of Hypertension with regard to the accuracy of the measured values.

Three measuring modes. Rapid single measurements, double measurements and manual measurements (auscultation with stethoscope).

Easy and quick hygienic reprocessing due to optimised design and without any moving control elements.

Electronic measurement technology exclusively manufactured for HEINE by ERKA.

EN200 BP Sets

EN200 BP as single module (long USB cord, E4-USB plug-in power supply), coiled tubing to connect cuffs, 3 Quick Cuff reusable cuffs (adult, small / adult / adult, large)	M-000.09.171
EN200 BP for use with the EN200 Wall Transformer (short USB cord), coiled tubing to connect cuffs, 3 Quick Cuff reusable cuffs (adult, small / adult / adult, large)	M-000.09.170

Accessories

Cuff Basket [01]	M-000.09.140
------------------	--------------



BLOOD PRESSURE CUFFS FOR EN 200 BP:

- **High quality metal quick connector.**
- Various cuffs for **every area of application.**
- **Excellent skin compatibility.** Latex and PVC free. Biocompatible.
- **Easy handling.** Cuff with integrated bladder. Simplified hygienic reprocessing.



Quick Cuff (Reusable)

Adult, small	20.5 – 28 Circumference/cm	M-000.09.803
Adult	27 – 35 Circumference/cm	M-000.09.804
Adult, large	34 – 43 Circumference/cm	M-000.09.805
Adult, XL	42 – 54 Circumference/cm	M-000.09.806

Quick Cuff Single-Patient (Single-Patient Use)

Adult, small (10 pcs.)	20.5 – 28 Circumference/cm	M-000.09.813
Adult (10 pcs.)	27 – 35 Circumference/cm	M-000.09.814
Adult, large (10 pcs.)	34 – 43 Circumference/cm	M-000.09.815

HEINE Optotechnik GmbH & Co. KG, Kientalstr. 7, 82211 Herrsching, Germany
Tel. +49 (0) 8152/38-0, Fax +49 (0) 8152/38-202, E-Mail: info@heine.com, www.heine.com



A-000.00.216 e 10/16.
We reserve the right to change specification without notice.
© 2016 HEINE Optotechnik GmbH & Co. KG. All rights reserved.



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

