



# X-Ray Viewer Single Slimline



ORDER CODE: SBXSL

## Product Description

- LCD display
- 2.4cm thickness
- 100,000 hours working time
- No globes required
- Lightweight
- High luminance

## Type and dimensions

- SBXSL Single Bay - Dimensions L 47 x W 2.4 x H 55cm



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995

