



healthtec

Supporting you

1300 362 923

www.healthtec.com.au

SPECIALTY

Gynae Examination Table

Designed in conjunction with a Gynaecologist this table provides sufficient and effective pelvic tilt for examining the cervix.

It also aids in the ease of manipulation of speculums.

The low minimum height with no lateral movement allows for easy patient access and transfer while the maximum height can quickly be set to suit the Practitioner's requirements reducing back strain and injury.



Features

- Electric Height Adjustment with footswitch
- X-Frame Lift action minimises pinch and crush points
- Heavy duty individually lockable castors
- Durable powder coated rust resistant frame
- High quality anti-microbial vinyl & CFC free foam
- 5 year warranty on frame & electrics
- Australian designed & manufactured

Specifications

- 250kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 450mm Minimum Height
- 950mm Maximum Height
- 710mm Table Width
- 1950mm Table Length
- Manual Backrest Adjustment +70°
- Electric pelvic seat section +20°



The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6 © October 2013



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995








Supporting you

PRICE SHEET

www.healthtec.com.au

Gynae Examination Table

| | CODE | QTY |
|--|------------------------|-----|
| Product Description | | |
| LynX Gynae Examination Table (710 wide) w/castors | 53751T-7 | |
| Width Options | | |
| 770mm Width Option | -75 | |
| Custom Width Option | 6014 | |
| Foam & Vinyl Options | | |
| 75mm Comfort Foam | 6005 | |
| Non-Stock Vinyl Option | 6018 | |
| Leather Upgrade | 6008 | |
| *please refer to our colour chart for vinyl colour options | | |
| Head Section Options | | |
| Nose Hole – Examination Backrest | 6011-2 | |
| Nose Hole Plug | 6010 | |
| Footswitch Option | | |
| Additional Footswitch (2 function) | 6024 | |
| Examination Options | | |
| Gynae Leg supports | 5114 | |
| Castor Options | | |
| Central Locking Castors | 6037 | |
| Hoist Access Clearance (140mm clearance) | 6028 | |
| Table Protection | | |
| Paper Towel Roll Holder | 6017 | |
| Disposable Fitted Sheet – 10 per pack | ANWFIT | |
| Extended Warranty | | |
| 10 years Frame & Electrics | 7800 | |
| Freight | | |
| Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options. | | |
| Shipping Dimensions | 2050mm x 850mm x 650mm | |
| Shipping Weight | 80kg (approximately) | |



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

