



Supporting you

1300 362 923

www.healthtec.com.au

TREATMENT TABLES

Two Section Treatment Table

Designed for the budget conscious practitioner looking for a reliable sturdy Australian made table built to last.

The manual adjustable head section can be slightly lowered allowing the prone patients neck muscles to relax giving optimum comfort.

Ideal for Physiotherapy, Massage or Osteopathy.



Features

- Electric Height Adjustment with footswitch
- X-Frame Lift action minimises pinch and crush points
- Available with either adjustable feet or castors
- Durable powder coated rust resistant frame
- High quality anti-microbial vinyl & CFC free foam
- 5 year warranty on frame & electrics
- Australian designed & manufactured

Specifications

- 250kg Maximum Lifting Capacity (moving)
- 400kg Maximum Working Load (stationary)
- 450mm Minimum Height
- 950mm Maximum Height
- 610mm Table Width
- 1840mm Table Length
- Manual Headrest adjustment +40° to -90°
- Nose Hole in headrest as standard

Available options:



5102: Procedure Arm



5105: Prone Arm Rest - Massage



6017: Paper Towel Roll Holder



The Healthtec Warranty

All Healthtec™ manufactured items are covered for defects in material or workmanship for 5 years from date of purchase. Original suppliers warranty applies.

TGA AR-00106892-6

© September 2013



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995








Supporting you

PRICE SHEET

www.healthtec.com.au

Two Section Treatment Table

Product Description	CODE	QTY
LynX Treatment Table – Two Section (610 wide) w/adjustable Feet	53020T	
LynX Treatment Table – Two Section (610 wide) w/castors	53021T	
Ezy-Move Wheel Option		
Retractable Ezy-Move Wheel	6000	
Width Options		
670mm Width Option	-65	
710mm Width Option	-70	
770mm Width Option	-75	
Custom Width Option	6014	
Foam & Vinyl Options		
75mm Comfort Foam	6005	
Non-Stock Vinyl Option	6018	
Leather Upgrade	6008	
*please refer to our colour chart for vinyl colour options		
Nose Hole Options		
Nose Hole – Extra (Footend)	6011	
Nose Hole Plug	6010	
Face Crest Pro (Footend)	5519	
Armrest & Footrest Options		
Procedure Arm	5102	
Prone Armrest (for facecrest pro)	5119	
Prone Armrest -Massage	5105-2	
Footrest Extension	5104	
Motor and Footswitch Options		
Electrically Operated Head Rest	6033	
Additional Footswitch (1 Function)	6020T	
Table Protection		
Paper Towel Roll Holder	6017	
Fitted Treatment Table Cover w/nose hole (Micro Fabric)	7201	
Extended Warranty		
10 years Frame & Electrics	7800	
Freight		
Freight charges are not included in our prices. Please contact your Healthtec representative for freight costs and options.		
Shipping Dimensions	2050mm x 850mm x 650mm	
Shipping Weight	80kg (approximately)	



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

