

V6

Vital Signs Monitor



- 🔗 New generation Vital Signs Monitor
- 🔗 7" high resolution display for easy reading
- 🔗 Lightweight, portable design and user-friendly interface for easy operation
- 🔗 Powerful storage capacity



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



General:

Size and Weight

- Size: 160×130×260mm (W×H×D)
- Weight: (2.5 kg) (Including a lithium battery)

Display

- Type: 7" Color Anti-glare TFT LCD
- Resolution: 234×480 pixels

Power supply

- Input voltage: AC: 100V~240V, 50Hz/60Hz
- Input power: 70VA
- Fuse: T1.6AL/250V, 2-Φ5×20mm

Battery

- Type: 11.1V 4000mAh rechargeable lithium ion battery
- Operating time: >8 hours
- Charge time: 6h to 100% (Standby)

Data Storage

- Clinic Mode:
 - Patient quantity: 1000
 - Recording number: 16000
 - Single patient clinic record number: 1~16000
- Monitor Mode:
 - Patient quantity: 1
 - Trend data: resolution 1Min: can store 96 hours
 resolution 5 Min: can store 480 hours
 resolution 10 Min: can store 960 hours
- Alarm events: 1000
- NIBP measurement record: 5000

Environmental Specifications

- Operating
 - Temperature: (5~40) °C
 - Humidity : ≤85% (noncondensing)
 - Atmospheric Pressure: (860~1060) hPa
- Transportation and Storage
 - Temperature: (-20~+55)
 - Humidity : ≤93% (noncondensing)
 - Atmospheric Pressure: (500~1060) hPa

PR

- Measurement range: 20~250 bpm
- Resolution: 1 bpm
- Accuracy: ±1% or ± 1 bpm, whichever is the greater

Performance:

BLT Digital SpO₂

- Measurement range (0~100) %
- Resolution: 1%
- Accuracy: ± 2% (70%~100%)
 ± 3% (35%~69%)
- unspecified (0~34%)

Nellcor SpO₂ (Option)

- Measurement range: (0~100) %
- Resolution: 1%
- Accuracy: ± 2% (70%~100%) (adult)
 ± 3% (70%~100%) (neonate)
 ± 2% (70%~100%) (low perfusion)
- unspecified (0~69)%

NIBP

- Technique: Automaticoscillometry
- Range: Adult: 10~270mmHg
 Child: 10~235mmHg
 Neonate: 10~135mmHg
- Accuracy: Static ±3 mm Hg
- Clinic: Average error: ±5 mmHg
- standard deviation: ≤8 mmHg
- Unit: mmHg.kPa
- Pulse rate range: 40-240 bpm

Microstream EtCO₂ (Option)


- Measurement way: Infrared spectrum
- Measurement mode: sidestream
- Sample flow rate: 50 ml/min ± 10 ml/min
- Measurement range: (0~19.7) %, (0~150 mmHg)
- Resolution: 0.1% or 1mmHg
- Unit: %, mmHg, kPa
- Accuracy: ±2 mmHg 0 mmHg to 40 mmHg,
 ±5% of reading 41 mmHg to 70 mmHg,
 ±8% of reading 71 mmHg to 100 mmHg,
 ±10% of reading 101 mmHg to 150 mmHg,
 (760mmHg, temperature is 25 °C)
 (when RR >80 rpm, all the range is ±12% of reading)

Standard Configuration

NIBP,BLT digital SPO₂,PR, Lithium Battery

Option

Micro Stream EtCO₂, Nellcor SpO₂, instant Infra-Red ear thermometer



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

