



# Medisafe Plus 1000



## Why Choose us

Medisafe Plus is the original premium brand of vaccine fridge in Australia. It is designed and built in Australia and has been protecting vaccines and medical products since 2010.

The Medisafe Plus Family of medical fridges are purpose-built for vaccine storage and are designed and built in Australia. All our medical fridges follow the "Strive for 5" National Vaccine Storage guidelines

A Medisafe Plus medical/vaccine fridge is feature packed, resulting in a vaccine storage solution that gives you peace of mind – knowing your vaccines are stored in the safest possible manner.

Medisafe Plus can tailor our solution to end user needs to allow for additional probes or temperature recording systems through the manufacturing process.

Medisafe Plus cloud monitoring is now a standard on all Medisafe Plus units, with 12 months free access to the portal. The cloud monitoring acts as an automated data logger which allows you to view and print your temperature chart with Min and Max and also data in real time on any web enabled smartphone, tablet, or computer from any location.

All units also the ability to have a BMS system plugged into the back of the unit for hospital remote monitoring.



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300633787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995





# Medisafe Plus 1000



## Medisafe Plus 1000

is one of the most feature-packed models of vaccine fridges available in Australia.

Medisafe Plus 1000 medical fridge is a purpose-built vaccine fridge designed and built in Australia in line with the National Vaccine Storage guidelines "Strive for 5".

Medisafe Plus vaccine fridges are feature packed, resulting in a vaccine storage solution that gives you peace of mind – knowing your vaccines are stored in the safest possible manner.

The Medisafe Plus 1000 is QCPP accredited and TGA approved.

The cloud monitoring acts as an automated data logger which allows you to view and print your temperature data in real time, from any location, on any web enabled smartphone, tablet, or computer.

All units also the ability to have a BMS system plugged into the back of the unit for hospital remote monitoring.

## Features

- 24 months Parts and Labour warranty
- Dedicated 24-hour technical support
- QCPP cold chain accredited
- TGA approved
- Min/Max Controller
- Temperature Controller
- Dual temperature probes for ambient and simulated vaccine temperatures
- AQ Box Cloud Monitoring – with SMS and email alerts of power or temperature failure
- Audio/Visual alarms
- Lock and Key
- Adjustable shelves
- Low Energy Consumption

## Optional Extras

- UPS Unit with Back up battery
- Preventative Maintenance program with extended warranty

<b>DIMENSIONS (HxWxD)</b>	<b>2200 x 1130 x 745</b>
<b>GROSS WEIGHT</b>	<b>210 Kg</b>
<b>GROSS CAPACITY</b>	<b>980 Lt.</b>
<b>TEMPERATURE RANGE</b>	<b>2° to 8°C</b>
<b>NUMBER OF SHELVES</b>	<b>10</b>
<b>POWER CONSUMPTION PER 24Hr</b>	<b>Tec = 3.57 kWh MEPS = 3.29 kWh</b>
<b>REFRIGERANT</b>	<b>R290a</b>

Local distributor:



Medical and Surgical Requisites  
 33 Fulcrum St, Richlands, QLD 4077  
 Phone : 07 3859 2900  
 Email: medical@medsurg.com.au