

VACCINE FRIDGES

FKG 311 MEDISAFE PLUS

Purpose built QCPP accredited vaccine fridge. 311 and 371 models are purpose built for larger storage space. Electronic temperature control and monitoring system ensures consistent temperatures. Data logger included allowing quick and easy download and recording of temperature history. Auto-dialler compatible. Lock and key included.



FKG 311 - MEDISAFE PLUS

FEATURES :

QCPP accredited vaccine fridges, Electronic temperature control & monitoring system, visual and audible auto resetting alarms, Min/Max temperature recorder, noise level of 55 db, low energy consumption, easy to use computer Data Logger, Australian design, 24 months warranty, power failure alarm, 311, 371 and 410 come standard with an integral lock and key, internal evaporator guard and door fan switch, Data logger has a visual alert, 20% larger blue digital display, power supply 10 amp, and dedicated national service and help line 1800 658 120.

DIMENSIONS & WEIGHT:

Height	Width	Depth	Gross Weight (kg)	Net Weight (kg)
1630	595	595	82	78

SPECIFICATIONS :

Gross Capacity litres	Net Capacity litres	Temperature Range °C	Power Consumption per 24h, kWh	Frequency HZ	Refrigerant
306	281	+2 to +8	3.69	50	R134a

UPS (UNINTERRUPTED POWER SUPPLY):

You can make your Medisafe Plus Range, UPS Protected, by purchasing the battery pack separate, that provides emergency power to your Medisafe blood fridges, when there is a power shutdown. UPS will give you up to 4 to 6 hours of back-up power depending upon the size and model of your fridge. **Remote monitoring device is also available.**

Dimensions: UPS: H-440mm, W-90mm, D-465mm
Battery Pack: H-440mm, W-132mm, D-465mm

VESRT 15 with Battery Pack



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995