

## VACCINE FRIDGES

### HB80 MEDISAFE PLUS

Purpose built QCPP accredited vaccine fridge. The smallest and most convenient unit in our range. Electronic temperature control and monitoring system ensures consistent temperatures. Data Logger included allowing quick and easy download and recording of temperature history. Auto-dialler compatible. Lock and key included.



HB80 - MEDISAFE PLUS

#### FEATURES :

QCPP accredited vaccine fridges, Electronic temperature control & monitoring system, visual and audible auto resetting alarms, Min/Max temperature recorder, noise level of 35.5 db, low energy consumption, easy to use Computer Data Logger, Australian design, 24 months warranty, power failure alarm, 311 and 371 come standard with an integral lock and key, Data Logger has a visual alert, 20% larger blue digital display, power supply 10 amp and dedicated national service and help line 1800 658 120

#### DIMENSIONS & WEIGHT:

Height	Width	Depth	Gross Weight (kg)	Net Weight (kg)
830	500	560	53	42

#### SPECIFICATIONS :

Gross Capacity litres	Net Capacity litres	Temperature Range °C	Power Consumption per 24h, kWh	Frequency HZ	Refrigerant
80	77	+2 to +8	3.08	50	R134a

#### UPS (UNINTERRUPTED POWER SUPPLY):

You can make your Medisafe Plus Range, UPS Protected, by purchasing the battery pack separate, that provides emergency power to your Medisafe blood fridges, when there is a power shutdown. UPS will give you up to 4 to 6 hours of back-up power depending upon the size and model of your fridge. **Remote monitoring device is also available.**

**Dimensions:** UPS: H-440mm, W-90mm, D-465mm  
Battery Pack: H-440mm, W-132mm, D-465mm

#### VESRT 15 with Battery Pack

