

PRODUCT SPECIFICATION SHEET

| | |
|---|--|
| Product Code: | MS-3-05 |
| Product Description: | Face Mask, Procedural, Ear-Loop, Blue, Type IIR (Level 2) |
| <u>Product Information</u> | |
| Type: | Procedural Face Mask, Type IIR (Level 2) |
| Material - Mask: | PP Non-Woven / PP Meltblown |
| Colour: | Blue |
| Mask Size: | 175 x 95mm |
| Design & Features: | Triple pleat design, 2x Ear-Loops, 3-Ply, Latex-Free, Non-Sterile |
| Barrier Protection Level: | EN 14683:2019 Type IIR (AS 4381: 2015 Level 2) |
| Bacterial Filtration Efficiency (BFE) %: | ≥ 98% |
| Differential Pressure (Delta P) : | <60 Pa/cm ² (<6.0 mm H ₂ O/cm ²) |
| Synthetic blood penetration resistance: | ≥ 16.0 kPa (>120mm Hg) |
| Intended Use: | Suitable for healthcare use where the risk of blood or body fluid splash is low to moderate. Single use only |
| Quality Assurance: | ARTG 339951 ISO 13485:2016, ISO 9001:2015 |
| Shelf Life: | 3 years from date of manufacture |
| Storage: | Store in a cool, dry area, away from direct sunlight and humidity |
| Country of Manufacture: | China |
| <u>Packing Information</u> | |
| Packing Description: | 50 pieces/box, 40 boxes/carton |
| Carton Quantity: | 2,000 pieces |

Created: 17/07/20
Last edited: 17/07/20



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

