


# Minispir®


mini-laboratory  
for spirometry and oximetry

# WinspiroPRO®

high performance  
PC software



**Pediatric Incentive Animations**



Plugs directly into the USB port.  
Real time Flow/Volume loop and Volume/time curve with PRE/POST comparison.  
Advanced spirometry test interpretation.  
Pediatric incentive animations.  
Estimated Lung Age (ELA).

Bronchial provocation test including new protocol with FEV1 response curve. **Mannitol**  
Temperature sensor for BTPS conversion.

**Option available: Oximeter module (it can be purchased separately)**

**WinspiroPRO®** is a unique spirometry and oximetry software, which comes standard with **Minispir®**

All patient records are shown on simple, single-screen patient cards with dynamic management of all data and graphs.

**WinspiroPRO®** can easily be connected to a hospital database or EMR and occupational health system. **(HL7 interface)**

**Supports NHANES III standard. Network Version available on request**



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

