

*One Step Ahead in Electrocardiography*

## Resting ECG




# 12 lead ECG

Comprehensive solutions for Mobile, Table Top or Hand Held ECG.

No computer experience is needed, due to the system's user friendly software, which allows even first-time users to retrieve, view and print ECG.

The system works on WINDOWS® environment or Android.

- Simple interface to local Hospital and Patient Information System
- Fax or E-mail ECG results directly from the PC or Android system
- Automatic Interpretation
- Caliper QRS and ECG measurements
- Comprehensive display and printout
- Laser or ink printout
- Physician remarks and signature
- Network abilities, data sharing
- Multiple languages
- Export formats: PDF, SCP, JPEG, XML, and GDT, HL7
- DICOM Interoperability




## Late Potentials - SAECG

- Detection of High Frequency WIDTH, RMS of last 40ms, duration under 40uV, High frequency Noise
- Complies with 1991 task force committee

## Heart Rate Variability

- Time and Frequency Domain analysis designed for Short Study.
- Measured/Reported parameters according to HRS/ESC Guidelines



[www.norav.com](http://www.norav.com) | [info@norav.com](mailto:info@norav.com)



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995




## ECG LAN (Green) - Mobile Wireless Network Solution

- Environmentally Friendly (paperless)
- Ideal for Hospitals
- Touch Screen
- Network Printing



## BLUE-ECG - Wireless Resting ECG

- Eliminates the bulky and expensive interface cable
- Caregivers have greater freedom of movement while acquiring ECGs.
- No cable induced artifacts
- Android, Smartphone or tablet compatible
- Windows compatible



## 1200M - The Classic model

- Ideal for physician offices
- PC-ECG connects via USB to any PC



### Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



## MEANS Interpretation

- World renowned automatic ECG Diagnosis –Adult and pediatric, developed at Erasmus University, the Netherlands.

## I200HR - High Resolution model

- Offers 16 bit resolution and very high sampling rate
- Ideal for hospitals, research laboratories, Signal Averaging and Heart Rate Variability applications



## NECG Thermal Printer Electrocardiographs

- Up to 4 hours battery life
- Automatic charging of the Lithium Ion battery when connected to main electricity
- Simultaneous 12 lead analysis with interpretation
- Continuous impedance test with alarm display
- Transfer of stored data to the Norav ECG Management System NEMS



NECG-3



NECG-12C  
12 Lead



NECG-12  
12 Lead

[www.norav.com](http://www.norav.com) | [info@norav.com](mailto:info@norav.com)



Established 1922


### Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 [www.medsurg.com.au](http://www.medsurg.com.au)  
Phone 1300medsurg (1300 633 787) Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au) Fax 07 3859 2995



## NEMS Cloud

- Upload ECGs from :  
Tablet/Smartphone or from a PC at remote clinic




## Norav ECG Management System NEMS:

Norav 12 leads ECG Stress system can work either as a stand alone application or it can be linked to Norav Medical ECG Management System (NEMS).

The NEMS manages studies and patients' data of:

- Resting ECG, Stress Test ECG, Long Term ECG Monitoring, Telemetry ECG, Heart Rate Variability, Late Potentials (SAECG), Holter ambulatory ECG, Spirometry Ambulatory Blood Pressure
- Easy interface to Hospital Information Systems (HIS) and Electronic Health Records (EHR)
- Comparison of studies



## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



See also our other product lines:

## Stress ECG



## Ambulatory ECG and Blood Pressure:



## Telemetry ECG

- 2,6 or 12 Lead Display
- Small Lightweight Transmitters
- Simple to Operate
- PC Based Central Station
- Arrhythmia Detection
- 12 Lead Surveillance Telemetry



CE 0473

[www.norav.com](http://www.norav.com) | [info@norav.com](mailto:info@norav.com)

290516