

NURSES KIT - NKSSAVC/006

Two-Hand Aneroid Sphygmomanometer with the Liberty Luxury Doctors Dual Head Stethoscope

Liberty Lightweight Aluminium Stethoscopes

Featuring a two-sided chest piece, only a truly extraordinary precision instrument combines outstanding performance with distinctive styling. Acoustic response is remarkable from its quality engineered chest piece. A very economical and well designed unit.

Stethoscope Features

- Light weight aluminium reversible head chest piece
- Adjustable chrome plated internal binaural
- Flexible 56 cm PVC tubing in 8 fashionable colours for departmental coding or personal identification
- Rubber non chill ring for patient comfort
- Stethoscope tubing length 76.2 cm
- Packaged in a hanging display case
- Name tag included

- Navy Blue
- Royal Blue
- Magenta
- Teal
- Purple



Colours are an indication only and may vary from actual product. Colours may vary between supplies.



Sphygmomanometer Features

- Black enamel 300 mmHg no stop pin manometer with attachment clip
- Velcro cuff, complete with bladder, bulb and valve
- Chrome plated deflation valve
- Two-handed use, suitable for both right and left hands
- Range 0-300 mmHg
- Complete with leatherette zipped case



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

