

Omron HEM907 Professional Blood Pressure Monitor



The Omron HEM-907 is the perfect addition to any clinical setting. With four unique operation modes, single mode for once off measurements, average mode for an average of 2 or 3 measurements, manual mode for auscultation method and check mode for comparison of pressure setting between the Omron HEM-907 and a sphygmomanometer. Clinically validated, the Omron HEM-907 is a safe, easy to use and accurate alternative to a sphygmomanometer.

- Clinically Validated
- 4 Unique Operation Modes
- 3 Cuff Sizes (S,M,L)
- Made in Japan
- Most widely used clinical BP monitor in Australian General Practices



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

