

OMRON

AUTOMATIC BLOOD PRESSURE MONITOR

HEM-7280T

Recognized with the iF Design Awards 2015



Wide Range Cuff
22 - 42cm

BRILLIANT DESIGN FEATURES WITH CONNECTIVITY



for iPhone / Android

For the list of supported smart devices
omronconnect.com/devices

For details, visit:
omronconnect.com

All for Healthcare



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



Established 1922

AUTOMATIC BLOOD PRESSURE MONITOR

HEM-7280T

Manage blood pressure with OMRON Automatic Blood Pressure Monitor HEM-7280T. Track your results with OMRON connect app.



Detects Irregular Heartbeat



Blood Pressure Level Indicator
For easy interpretation of results.



Average of last three readings with date and time



Morning Hypertension Indicator



Cuff Wrapping Guide
Indicates OK when cuff is wrapped snugly



Backlight
For easy reading



Wide Range Cuff
Fits arm circumference 22-42 cm



2 users X 100 memory



Intellisense Technology
Automatically applies the right amount of pressure for fast, accurate and more comfortable measurements.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



OMRON connect makes it easy to record, view and wirelessly sync your health data.



Simple to connect

Easy setup process for automatic data transfer



Easy Sync of Data

Share your health information with OMRON connect app



Clear And Insightful Results

Spot trends and monitor your health over time

For the list of supported smart devices
omronconnect.com/devices

For more details, visit:
omronconnect.com



TRACK YOUR HEALTH WITH EASE IN THE COMFORT OF YOUR HOME.

Many factors can influence your blood pressure readings; physical activity, time of the day, stress etc.

With the convenience of the OMRON Automatic Blood Pressure Monitor, you will be able to monitor your blood pressure regularly in the comfort of your home.

Trust your health to Omron Blood pressure Monitor, the Reliable partner for your Blood Pressure and health management needs!

*For people with high blood pressure. Consult your doctor to evaluate the readings.
Check your device periodically for accuracy.*



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



SIMPLE TO CONNECT

Easy setup process for automatic data transfer.

ACCESS ANYTIME

Your latest activity and personal history is always at your fingertips with your health data securely stored on your smartphone.

CLEAR AND INSIGHTFUL RESULTS


Spot trends and monitor your health over time.

INTUITIVE DASHBOARD


View recent measurements on the dashboard and get visual reminders to measure frequently.

KEEP TRACK OF YOUR HEALTH

See how your health has changed over time with clear graphs and charts.



Colour of screen lightens on Home Page: reminder for user to measure frequently



Graph format is correct at time of printing and subject to changes with upgrade of APP



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au


www.medsurg.com.au
Fax 07 3859 2995






OMRON connect

Wirelessly syncs measurements with your smartphone to easily view, save, and manage your health data.




Get Connected

Transfer measurements to your smartphone with ease after a simple, guided set-up.



Access Anytime

Your latest activity and personal health history are always at your fingertips through the intuitive dashboard.



See Progress

Health data is displayed in clear and insightful graphs that help you to see your health trends.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995



SPECIFICATIONS

Model	HEM-7280T
Measurement Method	Oscillometric Method
Measurement Range	Pressure: 0 to 299 mmHg Pulse: 40 to 180 beats/min
Accuracy	Pressure: ± 3 mmHg Pulse: $\pm 5\%$ of display reading
Battery Life	Approximately 1000 measurements (using new alkaline batteries)
Power Source	4 "AA" batteries 1.5V or Optional AC adapter
Console Weight	Monitor: Approximately 364g without batteries Arm Cuff: Approximately 170g
Cuff Circumference	22 to 42cm
Package Contents	Monitor, Arm Cuff, Instruction Manual, Storage Case, Battery Set

This OMRON product is produced under the strict quality system of OMRON HEALTHCARE Co., Ltd., Japan. The core component for OMRON blood pressure monitors, which is the Pressure Sensor, is produced in Japan.

Distributed by:

JA DAVEY PTY LTD
626 Lorimer Street Fishermans Bend,
VIC 3207, Australia

JA DAVEY LTD
C/- Healthcare Logistics 58 Richard Pearse
Drive Mangere, Auckland 2022,
New Zealand



OMRON HEALTHCARE Co., Ltd.
53, Kunotsubo, Terado-cho
Muko, KYOTO, 617-0002 JAPAN



OMRON HEALTHCARE EUROPE B.V
Scorpius 33, 2132 LR Hoofddorp
THE NETHERLANDS

Asia Pacific HQ

OMRON HEALTHCARE SINGAPORE PTE LTD.
438A Alexandra Road, #05-05/08
Alexandra Technopark, Singapore 119967
www.omronhealthcare-ap.com
CE0197



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

