

Planet ULED Examination Light with 6 x 1.6W LED Illumination

Evolving from the Planet U arm and 60 years experience of medical lighting manufacturing coupled with a decade of LED lighting applications, the ULED brings examination lighting into the future.

The average LED life (L70) is 50,000H+ which represents 15-20 years in normal use.

Features:

- available in 16,000 lux (ULED)
- quick release head
- IP64 light head shaped for easy cleaning/sterilisation
- total power use 12W
- low operational cost
- no maintenance required
- all metal construction for durability
- mounting brackets for every application (ceiling, wall and trolley mount)
- minimum intrusion into the workspace
- multi axis head rotation
- light weight - only 5kg
- designed and manufactured in Australia
- complies with AS/NZ 3200.241/IEC60601-241 (Medical electrical equipment - Particular requirements for safety - Surgical luminaires and luminaires for diagnosis)
- prompt supply of spare parts, technical service and support
- 5 year warranty

The ULED is suitable for:

- close examination
- over bed
- intensive care
- recovery rooms
- dermatology
- aged care



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

