

2 section couch



Name: 2 section couch

ORDER CODE:

ET2NB - Navy Blue

ET2BL - Black

ET2GR - Grey

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section through dual gas operated struts from -20o to +80o
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m2

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant

Specification:

- Total Length: 195cm
- Width: 71cm
- Head section length: 68cm
- Leg section length: 124cm
- Height: adjustable from 48cm to 90cm

Loading weight: 250Kgs



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

