

# 3 section couch medical



**Name:** 3 section couch medical

**ORDER CODE:**

**ET33NB - Navy Blue**

**ET33BL - Black**

- Electric Hi Lo operation by foot switch
- Sleek design with rounded frame legs for safety
- Minimum height 50cm to maximum height 90cm in 30 seconds
- Adjustable head section featuring double gas strut, angle adjustable from -45o to +70o
- Adjustable backrest featuring double gas struts, angle adjustable from -45o to +70o
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm castors for safety and mobility
- Retractable dual locking castors on all 4 corners
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m2

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant



## Specification:

- Total Length: 190cm
- Width: 71cm
- Head section length: 62cm
- Middle section length: 62cm
- Leg section length: 62cm
- Height: adjustable from 50cm to 90cm

Loading weight: 250Kgs



Established 1922

## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

