

Emergency Trolley



Name: Emergency Trolley

ORDER CODE:

EMT

Emergency Trolley

Dimensions:

650 x 480 x 980mm

Product Description:

- Defibrillator shelf
- Extension table
- 304 Stainless Steel IV Pole
- Cardiac CPR Board
- Waste Bin
- External Power Holder & Plug Board
- Oxygen Tank Holder
- 304 Stainless Steel Wire Basket
- ABS Divider
- Central or Seals Lock
- Castors With Brakes
- Non Slip Mat
- Luxurious Push Handle
- Soft Closing Drawers
- Drawer sizes: 3 x 110mm, 1 x 135mm, 1 x 195mm
- G.W 55kg
- Trolley has central locking system through side key

*All accessories included



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

