

Gynaecology table with leg extensions

Name: Gynaecology table with leg extensions

ORDER CODE:

GCNB - Navy Blue



- Operation by dual foot switch
- Electric Hi Lo with the other motor for adjustable seat angle
- Sleek design with rounded frame legs for safety
- Minimum height to maximum height in 30 seconds
- Fully adjustable / removable patient leg supports, also includes a plug for examination
- Gas struts control the back section
- Adjustable angle backrest from horizontal to +80°
- Adjustable angle seat from horizontal to +30°
- Excellent height range of 60cm - 100cm
- Standard width of 60cm
- Fitted with stainless debris tray
- 7.5cm castors for safety and mobility
- Retractable wheel system for ease of mobility
- T Motion motor
- Padding thickness 5cm
- Foam Density: 36g/m²
- CE approval

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant

Loading weight: 250Kgs



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

