

Podiatry / Multi Purpose Chair



Name: Podiatry / Multi Purpose Chair

ORDER CODE:

PCNB - Navy Blue

- Electric Hi Lo operation by foot switches
- Sleek design with rounded frame legs for safety
- Minimum height 55cm to maximum height 90cm in 30 seconds
- Dual motors. 1. Electric Hi Lo. 2. Electric mid-section
- Adjustable backrest featuring double gas struts, angle adjustable from 0o to +60o
- Adjustable seat section featuring one motor, angle adjustable from 0o to +20o
- Adjustable foot section featuring single gas strut
- Double removable arm supports
- Able to lift up to 250kgs working weight from its lowest height
- Retractable castors 50mm for excellent mobility
- Table goes horizontal for people recovery in emergency situations
- T Motion motors
- Padding thickness 5cm
- Foam Density: 36g/m²

Certification: CE, FDA, ISO9001-2008, ISO13485-2007

Upholstery: Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant

Specification:

- Backrest length: 90cm
- Seat section length: 56cm
- Leg section length: 37cm
- Foot section length: 20cm
- Upholstery width: 52 - 84cm
- Adjustable height: 55cm to 95cm
- Foam thickness: 5cm
- Foam density: 62cm
- Open: 220 x (52-84) x (55-95)cm

Loading weight: 250Kgs



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

