

Pacific Medical Ultrasonic Cleaners

Ultrasonic Cleaner 6 Litre



Product Code: PM-UCD06

- Turbo power: Normal and soft function to resolve the cleaning blind spots
- Auto-degas function for laboratory purposes
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel make it easy to operate
- Ultrasonic frequency: 40KHz
- All units come with basket and lid as standard

Capacity: 6L

Tank Size: 330 x 180 x 310 (LxWxH)

Heating power: 300W

Ultrasonic power: 150W

Temp Adjustable: ambient to 80°C

Drainage Valve: Yes

Digital Display

Timer: 1 - 99min



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Ultrasonic Cleaner 9 Litre



Product Code: PM-UCD09

- Turbo power: Normal and soft function to resolve the cleaning blind spots
- Auto-degas function for laboratory purposes
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel make it easy to operate
- Ultrasonic frequency: 40KHz
- All units come with basket and lid as standard

Capacity: 9L

Tank Size: 330 x 270 x 310 (LxWxH)

Heating power: 300W

Ultrasonic power: 200W

Temp Adjustable: ambient to 80° C

Drainage Valve: Yes

Digital Display

Timer: 1 - 99min



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Ultrasonic Cleaner 13 Litre



Product Code: PM-UCD13

- Turbo power: Normal and soft function to resolve the cleaning blind spots
- Auto-degas function for laboratory purposes
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel make it easy to operate
- Ultrasonic frequency: 40KHz
- All units come with basket and lid as standard

Capacity: 13L

Tank Size: 360 x 330 x 310 (LxWxH)

Heating power: 400W

Ultrasonic power: 300W

Temp Adjustable: ambient to 80°C

Drainage Valve: Yes

Digital Display

Timer: 1 - 99min



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Ultrasonic Cleaner 20 Litre



Product Code: PM-UCD20

- Turbo power: Normal and soft function to resolve the cleaning blind spots
- Auto-degas function for laboratory purposes
- Automatic ultrasonic frequency tracking
- Memory function
- LED display for temperature and timer
- User friendly and clear panel make it easy to operate
- Ultrasonic frequency: 40KHz
- All units come with basket and lid as standard

Capacity: 20L

Tank Size: 550 x 330 x 310 (LxWxH)

Heating power: 500W

Ultrasonic power: 400W

Temp Adjustable: ambient to 80°C

Drainage Valve: Yes

Digital Display

Timer: 1 - 99min



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

