

ELECTRIC PLASTER CUTTING SAW 125W INPUT POWER AND ALLOY BLADE

Liberty[®]
HEALTH PRODUCTS

new
design



Features

- Swing: Frequency: 133Hz
Range: ≥ 4 degrees
- Power Supply: 220V AC, 50Hz
- Input Power: 125W
- Blade Diameter: 62MM
- Power Cord Length: 2M
- Low noise total optimized mechanical structure and parts design <85dB
- Alloy blade is suitable for constant use
- Function of parenchyma protection is safe and reliable
- Variable speed control
- Suitable for cutting macromolecule plastic material

Complete With

- 1 Liberty[®] Electric Plaster Cutting Saw
- 2 Alloy Cutting Blades 62MM
- 1 Allen Key
- 1 Lockable Attache Style Carry Case

CODE	DESCRIPTION
PSEL	LIBERTY PLASTER SAW ELECTRIC
PSELB	LIBERTY PLASTER SAW SPARE BLADE



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

