

Innovative Medical Solutions

SUPPORTING THE SURGICAL & MEDICAL PROFESSION

Electronic Ear Irrigation

Announcing
the arrival of
the ProJet 101
in Australia



ProJet

The New Standard in Electronic Ear Irrigation



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995



Outstanding
value



Features:

- Power/ water pressure control button
- Water pressure indicator

Ear Irrigation kit consists of:

- Irrigator machine
- Water dispensing handle
- Water reservoir
- 10 ear irrigation tips
- Tip/ instrument storage
- Mains power cable and plug
- Foot switch

Specifications:

- Power: 230V, 50Hz
- Static pressure: 0.25 - 3 bar maximum
- Working time: 2 - 3 mins for 500cc
- Continuous usage: 8 mins
- Reservoir: 500ml
- Temp range: 5°C - 35°C
- Pressure: 860 - 1060 hpa
- HxWxD: 180x175x145mm

Utilising the most modern production methods the Guardian Proj 101 was designed to fulfil the most demanding surgeries and hospital ward requirements for a reliable, leak free and easy to use machine that is not expensive to run and does not require constant maintenance and subsequent down time.

The P101 has many unique features designed to help the busy surgery such as:

- 1 Extra long mains lead and handle hose so machine can be sited where it is most convenient
- 2 Integral wall mount frees up valuable work space by wall mounting the Proj[™] however easily removed for cleaning etc
- 3 Sliding flow regulator conveniently mounted in handle. Water flow can be temporarily stopped by sliding regulator to off position or flow regulated by sliding control to desired position
- 4 Ultra thin atraumatic tips help visualisation and direction of water flow
- 5 Digital control of water pressure and flow is easy to operate and remains constant and reproducible
- 6 Flat control panel is easy to keep clean
- 7 Tips held securely yet very easy to remove and replace
- 8 On board tip and/or instrument storage
- 9 Quiet operation
- 10 Oversize water reservoir silicone washer to reduce incidence of leaks
- 11 Economically priced

Water dispensing handle



Water reservoir



Proj[™]



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

