

Seca 213 Portable stadiometer

For measuring children and teenagers

The Seca stadiometer is particularly suitable for mobile use, e.g. for measuring children and teenagers in the course of medical examinations at schools. The measuring rod which can be dismantled into several pieces is easy to transport and can be set up anywhere in no time at all. A stable floor plate provides the necessary stability and the easy to read scale along both sides of the measuring rod means that the height can be read easily.

The Seca 213 is well-suited for medical examinations at schools, health centres, doctors surgeries and child health clinics

- Measurement range 20 - 207cm
- Graduation 1mm

Comfortable transport with measuring rod which can be dismantled into several pieces, as well as carrying handle and mounting option of the measuring rod parts on the floor plate.



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

