

Spirometry Mouthpieces

What mouthpiece suits what unit?

MicroGuard II Oval Bacterial Filter -White

SZ-CFN-V-892386

Suits:

- Carefusion Microlab, Microloop and SpiroUSB (with adaptor)
- QRS Orbit Disposable Mouthpieces
- MIR Disposable Turbine + Mouthpieces



Mada Bacterial Blue Filter

SZ-17530888

Suits:

- Carefusion Microlab, Microloop and SpiroUSB (without the use of a reusable adaptor)



Vitalograph Bacterial Filter

SZ-V28350

Suits:

- Vitalograph
- CosMed Pony and MicroQuark



Mada Bacterial Green Filter

SZ-175-30889

Suits:

- CosMed Pony and MicroQuark



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

