

Basic Task Lamps

- Even, diffused flood illumination
- Accepts incandescent or fluoro energy saver globes
- Globe not included
- Maximum reach 1070mm

Model I490 Superlux B Light

Unrestricted spring tension mechanism
Ideal for mounting on mobile base

Model I491 Superlux B Light with Wall Protection

Incorporates stop mechanism to prevent lamp damaging wall
Wall mount to be purchased separately



Mounting Options

Model I496	Wall Mount – Plastic moulded For models I491, I492, I493
Mount I497	Desk Clamp - Powdercoated For model I490, I492, I493
Mount I580	Mobile Base - Powdercoated For models I490, I492, I493

Magnifying Lamps

- 6 Diopter magnification through a 120mm diameter lens
- 10Watt circular LED produces shadowless 'daylight' illumination
- Angle poise spring balanced adjustable arm. Full extension 1200mm
- With protective dust cover flap
- Comes complete with desk clamp
- Mobile base or wall mount to be purchased separately

Model I493LED

Superlux Magnifying LED Lamp



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

