

# TD-1261D

## EAR INFRARED THERMOMETER

### FEATURES

- Fast ear scan in one second
- Dual-color backlight
- Probe cover ejection and detection
- Memory capacity up to 10 results



### THERMOMETER SPECIFICATIONS

|                               |                                                                                                                                            |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Temperature Measurement Sites | Ear Temperature                                                                                                                            |
| Temperature Range             | 32°C - 43°C                                                                                                                                |
| Accuracy                      | ±0.2°C (±0.4°F) for the measurement at 36°C-39°C (96.8°F-102.2°F)<br>±0.3°C (±0.5°F) outside above temperature range                       |
| Display Resolution            | 0.1°C / 0.1°F                                                                                                                              |
| Temperature Unit              | °C or °F set by factory                                                                                                                    |
| Backlight Screen              | <div> <div>32°C</div> <div>37.9°C</div> <div>43°C</div> </div> <div> <div>Lo</div> <div>Normal</div> <div>Fever</div> <div>Hi</div> </div> |
| Probe Cover Detector          | Equipped                                                                                                                                   |
| Memory Capacity               | 10                                                                                                                                         |
| Power Saving                  | 3 min                                                                                                                                      |
| Power Source                  | 2 x AAA                                                                                                                                    |
| Battery Life                  | 1000 measurements                                                                                                                          |
| Dimension                     | 164.8 (L) x 32.4 (W) x 38.7 (H) mm                                                                                                         |
| Data Transmission             | Bluetooth smart (TD-1261C); No Data Transmission (TD-1261D)                                                                                |
| Weight                        | 85.5g                                                                                                                                      |
| Units per Carton              | 20                                                                                                                                         |
| Net Weight per Carton         | 9.6 kg                                                                                                                                     |
| Gross Weight per Carton       | 10.6 kg                                                                                                                                    |
| Carton Size                   | 46.7 (L) x 36.7 (W) x 24.2 (H) cm                                                                                                          |

### PROBE COVER SPECIFICATIONS

|                            |                                   |
|----------------------------|-----------------------------------|
| Probe Cover Tray Dimension | 148.76 (L) x 46 (W) x 21.7 (H) mm |
| Probe Cover Unit per Box   | 320 pcs                           |
| Probe Cover Box per Carton | 6 pcs                             |
| Probe Cover Carton Size    | 31.5 (L) x 30 (W) x 15 (H) cm     |

**Mb<sup>®</sup>H**  
Mobile Health



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

