

# VanishPoint<sup>®</sup> Syringe

The passive safety mechanism of VanishPoint<sup>®</sup> minimises the risk of needlestick injuries



BEFORE

AFTER

- Activation cannot be by-passed
- Minimal training required
- Needle automatically retracts into the syringe

When safety matters,  
think TERUMO



Engineered for protection

**TERUMO**



**Medical and Surgical Requisites Pty Ltd**

Call within Brisbane 07 3859 2900  
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077  
Email [medical@medsurg.com.au](mailto:medical@medsurg.com.au)





[www.medsurg.com.au](http://www.medsurg.com.au)  
Fax 07 3859 2995




# VanishPoint® Syringe

LATEX FREE















## U-100 Insulin Syringes

Product Code	Colour	Volume	Gauge	Length		Shelf Pack	Carton
		mL	G	inches	mm	box	Qty
8RTI*15271		0.5	30	5/16	8	100	800
8RTI*15221		0.5	30	1/2	12.7	100	800
8RTI*10271		1	30	5/16	8	100	800
8RTI*10211		1	29	1/2	12.7	100	800

## Turberculin Syringes

Product Code	Colour	Volume	Gauge	Length		Shelf Pack	Carton
		mL	G	inches	mm	box	Qty
8RTI*10151		1	25	5/8	16	100	800
8RTI*10161		1	25	1	25	100	800
8RTI*10131		1	27	1/2	12.7	100	800

## Hypodermic Syringes

Product Code	Colour	Volume	Gauge	Length		Shelf Pack	Carton
		mL	G	inches	mm	box	Qty
8RTI*10371		3	20	1	25	100	600
8RTI*10381		3	20	1 1/2	38	100	600
8RTI*10361		3	21	1 1/2	38	100	600
8RTI*10311		3	23	1	25	100	600
8RTI*10321		3	23	1 1/2	38	100	600
8RTI*10301		3	25	5/8	16	100	600
8RTI*10391		3	25	1	25	100	600
8RTI*13011		3	25	1 1/2	38	100	600
8RTI*10571		5	20	1	25	100	600
8RTI*10581		5	20	1 1/2	38	100	600
8RTI*10551		5	21	1	25	100	600
8RTI*10531		5	22	1	25	100	600
8RTI*10541		5	22	1 1/2	38	100	600
8RTI*11041		10	22	1 1/2	38	100	600

VanishPoint® is a Registered Trademark of Retractable Syringe Technologies Inc. GHP140701D659V2 | PBD00135

### Terumo Australia Pty Ltd

Macquarie Park NSW 2113 Australia  
 Mount Albert Auckland 1025 New Zealand  
 E: cs\_australia@terumo.co.jp  
 Orders: au\_orders@terumo.co.jp

T: 1800 837 866  
 T: 0800 66 77 57



## Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900 33 Fulcrum Street Richlands Queensland 4077 www.medsurg.com.au  
 Phone 1300medsurg (1300 633 787) Email medical@medsurg.com.au Fax 07 3859 2995

