

WOLFMED Finger Pulse Oximeter



DZ-WM310

DESCRIPTION

- Support 6 display modes.
- Dual colour OLED display.
- Lightweight and Easy-to-use.
- Visual and sound reminder function.
- Algorithm inside.
- Low perfusion: 0.3%.
- Anti-Movement.
- Adjust the parameters in friendly menu.
- Ip22 water splash proof.
- Low battery voltage indicator.
- Automatically switch off while no signal.
- Standard two AAA 1.5V alkaline batteries support more than 20 hours continuous work.



Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

