

ProShield®

Mask Range

Fluid
Resistant

PROSHIELD® N-95

PARTICULATE RESPIRATOR MASK

- Protection against TB
- NIOSH approved
- Fluid resistant (160mmHg)
- AS/NZ 1716:2012 Respiratory protective devices

THERAPIES. HAND IN HAND.

BSNmedical


an Essity company

PROSHIELD® N-95

- Meets NIOSH criteria for protection against TB
- Secure full facial fit
- Easy donning
- Conformable nose piece
- Durable and comfortable headbands
- Optimum protection against blood and body fluid splash (160mmHg)



Features and benefits



BFE

With a **Bacterial Filtration Efficiency (BFE)** of greater than 99% for media greater than 3 microns this mask provides superior protection from potential infection.¹

PFE

The super high **Particulate Filtration Efficiency (PFE)** material filters >99% of particles greater than 0.1 microns. When securely fitted, the mask helps to provide an effective barrier to particulate matter from laser plumes and electrocautery procedures.²

N

The N95 mask is NIOSH approved as a N95 particulate filter respirator. It meets or exceeds the standard performance criteria demanded by NIOSH for the management of Tuberculosis (TB).³



The Fluid Resistant qualities of the mask provide protection against fluid strikethrough.⁴

File Notes:

1. BFE (Bacterial Filtration Efficiency %) is an independent test conducted by Nelson Laboratories that measures the effectiveness of face mask materials to filter out bacteria.
2. PFE (Particulate Filtration Efficiency) measures the effectiveness of the mask material to filter particulate matter.
3. NIOSH - National Institute of Occupational Safety and Health (USA).
4. Fluid resistance is a standard test for resistance of medical face masks to penetration by synthetic blood.

ORDERING INFORMATION

Ref-no.	Description	Items per unit	Shipper Qty
72509-09	ProShield® N-95 Small	bx/50	6 bxs
72509-10	ProShield® N-95 Medium	bx/50	6 bxs



Established 1922

Medical and Surgical Requisites Pty Ltd

Call within Brisbane 07 3859 2900
Phone 1300medsurg (1300 633 787)

33 Fulcrum Street Richlands Queensland 4077
Email medical@medsurg.com.au

www.medsurg.com.au
Fax 07 3859 2995

